



Name: _____ Employee ID #: _____

Department: _____ Site: _____ Date: _____

True or False (use “T” for true and “F” for false):

- _____ 1. The presence of excess amounts of blood or saliva in the collected specimen can interfere with test results.
- _____ 2. The performance of the internal procedural controls is documented with each patient and quality control result.
- _____ 3. Negative results can occur from inadequate specimen collection or antigen level below the detection limit of the test.
- _____ 4. An Extraction Reagent bottle can be used with a Test Device from a different kit lot number.
- _____ 5. Test specimens heavily colonized with *Staphylococcus aureus* can yield false positive results.
- _____ 6. The ICON® SC Strep A Rapid Test can detect non-viable and viable Strep A organisms; therefore the test may detect organisms which cannot be demonstrated in culture.
- _____ 7. The ICON® SC Strep A Test Device should be stored at 2-30° C in its original sealed pouch, out of direct sunlight.
- _____ 8. The test results are read after 10 minutes, but not after 15 minutes.
- _____ 9. Three (3) drops of Reagent A and three (3) drops of Reagent B are added to Extraction Well of the Test Device.
- _____ 10. All negative ICON® SC Rapid Strep A results are to be confirmed by culture.
- _____ 11. The use of external positive and negative controls assure that the test reagents are working properly and the user has performed the test correctly.
- _____ 12. Patient results can be reported when the internal procedural controls have not worked properly or when the external quality controls have failed.
- _____ 13. The date when external quality control is successfully performed is recorded on the JHMI QC Label affixed to the test kit box.

**The Johns Hopkins Medical Institutions
The Johns Hopkins Hospital Point-of-Care Testing Program
Rapid Strep A Testing (ICON® SC Strep A)
Test for Operators**

NAME: _____

Employee ID#: _____

Multiple Choice (circle all answers that apply)

1. Specimens for the ICON® SC Strep A test should be collected with:
 - a. A calcium alginate swab
 - b. A cotton-tipped swab
 - c. A rayon swab supplied with the kit
 - d. A Culturette swab with modified Stuart's transport media
 - e. All of the above

2. The operator is responsible for:
 - a. Proper storage of test kit and controls
 - b. Proper documentation of QC and remedial actions for failed QC.
 - c. Proper specimen collection
 - d. Adherence to test procedure to achieve optimal test results.
 - e. Documentation of internal procedural control performance with each patient and QC result.
 - f. All of the above

3. Specimens for the ICON® SC Strep A test are collected from:
 - a. The tonsils
 - b. The tongue and cheek surfaces
 - c. The back of the throat
 - d. Saliva from the gums and teeth
 - e. All of the above

4. Two levels of external quality control are performed on this test kit:
 - a. Monthly
 - b. When a new test kit box is opened
 - c. Weekly on each opened box
 - d. All of the above
 - e. None of the above

5. The internal procedural controls for the ICON® SC Strep A test include:
 - a. A reddish-purple Test Line that appears
 - b. The Extraction Reagent changes color
 - c. The background clears on the Test Device
 - d. A reddish-purple Control Line that appears
 - e. All of the above

6. Which of the following is **NOT** a correct interpretation of the Rapid Strep A test results?
 - a. Negative result – a reddish-purple Control Line but no reddish-purple Test Line
 - b. Invalid result – no reddish-purple Control Line and a reddish-purple Test Line
 - c. Positive result – a reddish-purple Control Line and two reddish-purple Test Lines
 - d. Positive result – a reddish-purple Control Line and a reddish-purple Test Line
 - e. Invalid result – a reddish-purple Control Line and dense background color