



## GeneXpert® Xpress System Maintenance Log

Analyzer ID(s): \_\_\_\_\_

Testing Unit: \_\_\_\_\_

Month: \_\_\_\_\_

Year: \_\_\_\_\_

**Refer to Section XV in the procedure for step-by-step instructions. At a minimum, each section must be initialed upon completion. Date all required areas as indicated.**

<b>Daily Maintenance (Decontaminate Testing Surface)</b>																														
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
<b>Weekly Maintenance (Power Down System to Clear Unwanted Temporary Files)</b>																														
Task	Week 1 <i>(Date and Initial Indicates Both Instruments Done)</i>					Week 2 <i>(Date and Initial Indicates Both Instruments Done)</i>					Week 3 <i>(Date and Initial Indicates Both Instruments Done)</i>					Week 4 <i>(Date and Initial Indicates Both Instruments Done)</i>					Week 5 <i>(Date and Initial Indicates Both Instruments Done)</i>									
Exit out of Xpress Software																														
Turn off Hub Computer																														
Power Down Hub																														
Power Down Instrument																														
Power Up Instrument and Hub																														
Login Cepheid-Admin (Password cphd)																														
<b>Quarterly Maintenance (Clean GeneXpert Xpress Instrument and Hub Surfaces)</b>																														
Task	Due Date: January 31 <i>(Date and Initial Indicates Both Instruments Done)</i>					Due Date: April 30 <i>(Date and Initial Indicates Both Instruments Done)</i>					Due Date: July 31 <i>(Date and Initial Indicates Both Instruments Done)</i>					Due Date: October 31 <i>(Date and Initial Indicates Both Instruments Done)</i>					As Needed <i>(Date and Initial Indicates Both Instruments Done)</i>									
Complete Tasks 1-4 of Weekly Maintenance																														
Use Ethanol Solution to Clean Exterior																														
Complete Tasks 5-6 of Weekly Maintenance																														