

## POC1-012F7a: TEMPERATURE LOG

(OraQuick® ADVANCE Rapid HIV-1/2 Antibody Test)

**Test Kit Storage Room**  
 Location: \_\_\_\_\_

**Control Kit Storage (Refrigerator)**  
 Location: \_\_\_\_\_  
 Refrigerator Serial #: \_\_\_\_\_

**Testing Room**  
 Location: \_\_\_\_\_

**Acceptable Temperature Range:**

**Test Kit Storage:** 2°C to 27°C  
 35°F to 80°F

**Control Kit Storage:** 2°C to 8°C  
 35°F to 46°F

**Testing Room:** 15°C to 37°C  
 59°F to 99°F

**Month/Year:** \_\_\_\_\_

Day	Temp. (°C) MIN	Temp. (°C) MAX.	Temp. (°C) CURRENT	COMMENTS / ACTIONS TAKEN	Initials	Reviewed By:
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						
11						
12						
13						
14						
15						
16						
17						
18						
19						
20						
21						
22						
23						
24						
25						
26						
27						
28						
29						
30						
31						

**Comment/Action Codes:**

- |   |  |
|---|--|
| 1. Inventoried/stocked/removed supplies               | 4. Ensured Room air intake/supply ducts not blocked      |
| 2. Removed overstock; ensured refrig. fan not blocked | 5. Called Maintenance- supplies/testing moved to: _____. |
| 3. Adjusted thermostat                                | 6. Continue to monitor                                   |

**Revision: 2/2011**