



POCT Coagulation Room Temperature Log

Year: _____

Unit/Room: _____

Month: _____

| DAY | Minimum ≥15°C | Maximum ≤30°C | Current Temp. (15°-30°C) | Initials | Comments/Actions Taken | Unit Review | Lab Review |
|-----|------------------|------------------|--------------------------------|----------|------------------------|----------------|---------------|
| 1 | | | | | | | |
| 2 | | | | | | | |
| 3 | | | | | | | |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | | | | | | | |
| 7 | | | | | | | |
| 8 | | | | | | | |
| 9 | | | | | | | |
| 10 | | | | | | | |
| 11 | | | | | | | |
| 12 | | | | | | | |
| 13 | | | | | | | |
| 14 | | | | | | | |
| 15 | | | | | | | |
| 16 | | | | | | | |
| 17 | | | | | | | |
| 18 | | | | | | | |
| 19 | | | | | | | |
| 20 | | | | | | | |
| 21 | | | | | | | |
| 22 | | | | | | | |
| 23 | | | | | | | |
| 24 | | | | | | | |
| 25 | | | | | | | |
| 26 | | | | | | | |
| 27 | | | | | | | |
| 28 | | | | | | | |
| 29 | | | | | | | |
| 30 | | | | | | | |
| 31 | | | | | | | |

NOTE: If outside of Temp. range, note actions taken in comment section and contact POCT Office. (5-2645)

Revision: 11/22/05

Comment/Action Codes:

1. Ensured Room air intake/supply ducts not blocked
2. Adjusted thermostat
3. Continue to monitor
4. Called Facilities/Maintenance- supplies moved to: _____ (room location)